

# FAQ Device BAB00

## 1. Terms of use

### 1.1 Calibration

**Is calibration necessary to use the Device BAB00?**

No, the device calibrates itself at predefined intervals. There are no manual steps involved in this process.

### 1.2 Usage environment

**Is laboratory environment and laboratory personnel required to use the Device BAB00?**

No, employees can operate the device on their own after a simple briefing. The Device BAB00 can be used in any enclosed space.

### 1.3 Service intervals

**Does the Device BAB00 require maintenance?**

No regular maintenance is compulsory, but we recommend a service interval of once a year.

## 2. Data storage

### 2.1 Sample data

**How much sample data can be saved?**

The system and software have no limitations. The number of data records that can be saved is only limited by the hard disk capacity of the computer used.

### 2.2 LIMS

**Can I connect my LIMS?**

Please feel free to contact us with your specific requirements for connecting your LIMS.